

Goretastic Fake Blood - Instructions

(40g sachet makes approx 2 pints of end product)

What you will need:

- A measuring jug with measurements in "ml"
- A kettle
- A large spoon
- A large aperture funnel
- An empty bottle or container to pour the blood into after mixing.

Instructions:

Add 200ml of cold water to a measuring jug, add in 600ml of boiling water to the jug and give a quick stir. Open the silver sachet and add in the powdered contents to the jug, stir vigorously until no powder is visible, keep stirring for 2 – 3 minutes until it thickens and to avoid a lumpy consistency.

Stir in additional cold water to make up to 1150ml of blood or until the consistency you require is met.

Allow to cool for approx 15 minutes. You're now ready for bottling!

Use your funnel to pour the mixture into the container / bottle.

Remember to wash your jug, spoon and funnel immediately after.

ONCE BLOOD IS MIXED AND COOLED USE WITHIN 2 DAYS, OR 5 DAYS IF KEPT REFRIDGERATED.

SAFE IN MOUTH BUT TRY TO AVOID INGESTION

Caution:

Children require adult supervision. May stain clothing and fabric. Will stain skin temporarily. Remove with soap and warm water. AVOID CONTACT WITH EYES.

ORDER MORE AT: <https://www.goretastic.co.uk>

Goretastic Horror Products is an Ace Group brand
Ace Group, 48 Brittain Avenue, Newcastle under Lyme, Staffordshire, ST5 7NT

Goretastic Fake Blood - Instructions

(20g sachet makes approx 1 pint of end product)

What you will need:

- A measuring jug with measurements in "ml"
- A kettle
- A large spoon
- A large aperture funnel
- An empty bottle or container to pour the blood into after mixing.

Instructions:

Add 100ml of cold water to a measuring jug, add in 300ml of boiling water to the jug and give a quick stir. Open the silver sachet and add in the powdered contents to the jug, stir vigorously until no powder is visible, keep stirring for 2 – 3 minutes until it thickens and to avoid a lumpy consistency.

Stir in additional cold water until the consistency you require is met or to create more liquid.

Allow to cool for approx 15 minutes. You're now ready for bottling!
Use your funnel to pour the mixture into the container / bottle.

Remember to wash your jug, spoon and funnel immediately after.

ONCE BLOOD IS MIXED AND COOLED USE WITHIN 2 DAYS, OR 5 DAYS IF KEPT REFRIDGERATED.

SAFE IN MOUTH BUT TRY TO AVOID INGESTION

Caution:

Children require adult supervision. May stain clothing and fabric. Will stain skin temporarily. Remove with soap and warm water. AVOID CONTACT WITH EYES.

ORDER MORE AT: <https://www.goretastic.co.uk>

Goretastic Horror Products is an Ace Group brand
Ace Group, 48 Brittain Avenue, Newcastle under Lyme, Staffordshire, ST5 7NT

Goretastic Fake Blood - Instructions

(10g sachet makes approx half pint of end product)

What you will need:

- A measuring jug with measurements in "ml"
- A kettle
- A large spoon
- A large aperture funnel
- An empty bottle or container to pour the blood into after mixing.

Instructions:

Add 50ml of cold water to a measuring jug, add in 150ml of boiling water to the jug and give a quick stir. Open the silver sachet and add in the powdered contents to the jug, stir vigorously until no powder is visible, keep stirring for 2 – 3 minutes until it thickens and to avoid a lumpy consistency.

Stir in additional cold water until the consistency you require is met or to create more liquid.

Allow to cool for approx 15 minutes. You're now ready for bottling!
Use your funnel to pour the mixture into the container / bottle.

Remember to wash your jug, spoon and funnel immediately after.

ONCE BLOOD IS MIXED AND COOLED USE WITHIN 2 DAYS, OR 5 DAYS IF KEPT REFRIDGERATED.

SAFE IN MOUTH BUT TRY TO AVOID INGESTION

Caution:

Children require adult supervision. May stain clothing and fabric. Will stain skin temporarily. Remove with soap and warm water. AVOID CONTACT WITH EYES.

ORDER MORE AT: <https://www.goretastic.co.uk>

Goretastic Horror Products is an Ace Group brand
Ace Group, 48 Brittain Avenue, Newcastle under Lyme, Staffordshire, ST5 7NT